


Purification and Detection for DYKDDDDK Tagged Proteins



GenScript's anti-DYKDDDYK resins offer high capacity, specificity, sensitivity, and reusability, making it ideal for affinity purification and immunoprecipitation of DYKDDDDK-tagged fusion proteins expressed in common systems such as bacteria, yeast, and mammalian cells.

High Purity

High Sensitivity

High Binding Capacity

Purity >95% under different expression system

Capable of binding as low as 5 ng of protein

>1.0 mg/ml resin


Anti-DYKDDDDK Affinity Resins Selection Guide

		Anti-DYKDDDDK Affinity Resin Easy NEW	Anti-DYKDDDDK G1 Affinity Resin	MonoRab™ Anti-DYKDDDDK Affinity Resin
Cat. No.		L00907	L00432	L00766
Volumes	Volumes	1ml, 5ml, 10ml, 25ml	1ml, 5ml, 10ml, 25ml	1ml, 5ml, 25ml
	Immunoprecipitation (IP)	√	N/A	ChIP
	Purification	√	√	√
	Binding capacity	~1.0 mg/ml resin	~1.5 mg/ml resin	~1.5 mg/ml resin
Elution	Alkaline elution	√	√	√
	Acid elution	N/A	\checkmark	N/A
	Neutral elution	N/A	\checkmark	N/A
	Peptide competition elution	$\sqrt{}$	\checkmark	√
	SDS-PAGE loading buffer	$\sqrt{}$	$\sqrt{}$	\checkmark
Regeneration times		4-10 times	4-10 times	4-10 times
Advantages		IP application: Low non-specific adsorption, High sensitivity Purification application: High purity, especially in yeast system	Purification application: High binding capacity, High purity	Tolerant to high salt and acid

Application Data for IP


Non-specific adsorption of HA tag fusion protein 1 2 3 4 Lanes: 1: Competitor S 2: Competitor T 3: Competitor ST 4: L00907 IP-WB pattern of non-specific adsorption of cGFP tag fusion protein 1 2 3 4 Lanes: 1: Competitor S 2: Competitor T 3: Competitor ST 4: L00907 IP-WB pattern of non-specific adsorption of HA and cGFP tag fusion protein

Conclusion: GenScript L00907 has almost no non-specific adsorption on HA and cGFP tag fusion protein, and its performance is better than that of competitors.




Conclusion: GenScript L00907 can still effectively capture the target protein under the experimental conditions with a very small amount of target protein, and its sensitivity is higher than that of competitors.

Application Data for Purification



Conclusion: The purity of the target protein purified by GenScript L00907 in different expression systems was >95%, especially in the *Yeast* system.



Anti-DYKDDDDK Detection Antibodies and Plates

Produc	ct Categories	Cat. No	Product Name	
	Unconjugated antibody	A00187	THE™ DYKDDDDK Tag Antibody, mAb, Mouse	
		A01428	THE™ DYKDDDDK Tag Antibody [HRP], mAb, Mouse	
Mouse monoclonal		A01429	THE™ DYKDDDDK Tag Antibody [Biotin], mAb, Mouse	
antibody:	Conjugated antibody	A01632	THE™ DYKDDDDK Tag Antibody [FITC], mAb, Mouse	
High sensitivity,	Conjugated antibody	A01809	THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse	
high specificity		A01810	THE™ DYKDDDDK Tag Antibody [iFluor 555], mAb, Mouse	
		A01811	THE™ DYKDDDDK Tag Antibody [iFluor 647], mAb, Mouse	
Rabbit	Unconjugated antibody	A01868	MonoRab™ DYKDDDDK Tag Antibody, mAb, Rabbit	
monoclonal		A01869	MonoRab™ DYKDDDDK Tag Antibody [HRP], mAb, Rabbit	
antibody:	Conjugated antibody	A01870	MonoRab™ DYKDDDDK Tag Antibody [Biotin], mAb, Rabbit	
High affinity		A01871	MonoRab™ DYKDDDDK Tag Antibody [FITC], mAb, Rabbit	
	Black	L00455B	DYKDDDDK Tag Antibody Plate (Black, 96-well)	
Detection Plate	Clear	L00455C	DYKDDDDK Tag Antibody Plate (Clear, 8X12 strip)	
	White	L00455W	DYKDDDDK Tag Antibody Plate (White, 96-well)	

